

# Resectr™ Tissue Resection Device

## Diagnose and treat, in one office visit

Remove multiple polyps under complete visualization with one entry into the uterus.

Resectr is the only oscillating resection device that works with your existing fluid management system, wall suction and hysteroscopes.

- ✔ Ready-to-use
- ✔ No complex set up
- ✔ Short learning curve
- ✔ Completely disposable
- ✔ Works with your existing equipment



### Low-tech, high returns

Streamline and accelerate the diagnosis and treatment of polyps, one of the most common causes of AUB, while maximizing reimbursement. Contact your Minerva Surgical sales representative for more information.

# Choose Resectr from Minerva for control without added complexity or costs.

## Compatible with appropriately-sized hysteroscopes




5Fr



9Fr

Each resecting blade is bi-directional, internally rotating and oscillating to provide six rotations per handle squeeze.

## See and treat multiple small polyps



Before



After

## Control at your fingertips

Squeeze and release the handle with your fingers to engage the resection device. Manual control means you can perform tissue resection based on what you see and feel during the procedure.



## Ordering information

Product Code	Description	Quantity
M0065907051	Resectr™ Tissue Resection Device 1.65mm / 5 French OD, 5mm window, and 36cm Long Cannula	Single
M0065907041	Resectr™ Tissue Resection Device 3mm / 9 French OD, 7.5mm window, and 36cm Long Cannula	Single



Caution: US Federal law restricts this device to sale by or on the order of a physician. The physician using the device must be trained in diagnostic hysteroscopy. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

All trademarks are the property of their respective owners.

Minerva Surgical  
4255 Burton Dr  
Santa Clara, CA 95054  
855-646-7874

[www.minervasurgical.com](http://www.minervasurgical.com)

© 2020 Minerva Surgical or its affiliates. All rights reserved.

K0084 Rev. A