

PRESS RELEASE

### **Minerva Surgical Introduces Fluid Deficit Readout Product for Operative Hysteroscopy**

Santa Clara, California, August 20, 2024, Minerva Surgical, a manufacturer and distributor of minimally invasive gynecologic surgical technologies, today announced the launch of the Symphion Fluid Deficit Readout, an optional accessory that provides fluid deficit volume readings during operative hysteroscopy procedures.

Operative hysteroscopy is a surgical procedure that requires precise fluid management to prevent complications. The Symphion Fluid Deficit Readout, which is intended for use with the Symphion Operative Hysteroscopy System, automatically calculates the fluid deficit during the procedure by measuring the weight of the saline bag throughout the procedure in real time and converting the weight measurement into fluid volume readings.

Additional key features of the Symphion Fluid Deficit Readout includes a fluid deficit accuracy within +/- 50mL, a user-friendly interface with an intuitive display for easy operation, and a compact design that allows for quick set up with the Symphion controller

“Continuous innovation with a focus on enhancing patient safety has always been the primary objective for Minerva Surgical. The Symphion Operative Hysteroscopy System is intended to volumetrically limit a patient’s exposure to fluid so the risk of fluid absorption and overload can be avoided. With the introduction of the Symphion Fluid Deficit Readout, the system has a new layer of procedural safety to complement the unprecedented precision during uterine cavity tissue resection.”, states Minerva Surgical Chief Medical Officer Eugene Skalny, M.D.

The Symphion Fluid Deficit Readout is the latest product launch to the Minerva Surgical portfolio of minimally invasive technologies. The company is committed to advancing gynecologic surgery and providing solutions that enhance outcomes.

**About Minerva Surgical, Inc.** Minerva Surgical is a commercial-stage medical technology company focused on developing, manufacturing, distributing, and commercializing minimally invasive solutions to meet the distinct uterine healthcare needs of women. The Company has established a broad product line of commercially available, minimally invasive alternatives to hysterectomy, which are designed to address the most common causes of Abnormal Uterine Bleeding (AUB) in most uterine anatomies. The Minerva Surgical solutions can be used in a variety of medical treatment settings and aim to address the drawbacks associated with alternative treatment methods and to preserve the uterus by avoiding unnecessary hysterectomies.

For more information about the innovative medical devices of Minerva Surgical, please visit [www.MinervaSurgical.com](http://www.MinervaSurgical.com).

**Contact:**

Minerva Surgical  
4255 Burton Drive  
Santa Clara, CA 95054  
Media/Press: [kevin.tracey@minervasurgical.com](mailto:kevin.tracey@minervasurgical.com)