

HERizon™

Detect, Diagnose, and Treat Uterine Cavity Conditions

16 degree cannula tip
with 100 degree field of view,
two LED lights in the tip

Thin 4.8mm
diameter sheath


2mm operative channel
for instruments up to 6 FR

Light 180g body
with ergonomic
design

Versatile monitor compatibility



Image Processor for video,
light and image capture




No cervical dilation

Minimal anesthesia

Reduce cross contamination

Compact office foot print

Hysteroscopy Accessories - Scissors, Forceps and Graspers



Scissors



Alligator



Foreign Body Grasper



Biopsy

To order contact Minerva Surgical Customer Service

By phone: call **855.646.7874**, and press for option 1, or email: customerservice@minervasurgical.com, or fax 866.465.2875

Product Name	Catalog Number	Quantity	Description
Medical Endoscope Image Processor	JY-MIP-1000	1 Unit	Camera and light source control unit
Disposable Electronic Hysteroscope	JY-DEH-001	1 Unit	Disposable 4.8 mm hysteroscope w/operative port, orders sold in multiples of 10 units
Disposable Endoscopic Forceps - Foreign Body Grasper	JY-C-350-18-S	1 Unit	For IUD grasping and removal under hysteroscopic guidance, orders sold in multiples of 10 units
Disposable Endoscopic Forceps - Alligator	JY-CC-350-18-S	1 Unit	For tissue biopsy or sampling under hysteroscopic guidance, orders sold in multiples of 10 units
Disposable Endoscopic Scissors	JY-JD-17-350-B	1 Unit	For cutting tissue under hysteroscopic guidance, orders sold in multiples of 10 units
Disposable Endoscopic Forceps - Biopsy	JY-PK-350-18-S	1 Unit	For tissue biopsy or sampling under hysteroscopic guidance, orders sold in multiples of 10 units

Intended use: The Disposable Electronic Hysteroscope is intended to be used for viewing of the adult cervical canal and uterine cavity for the purpose of performing diagnostic and operative procedures.

Contraindications: This product is not suitable for people with acute pelvic inflammatory disease. Hysteroscopy may be contraindicated by the following conditions, depending on their severity or extent: inability to distend the uterus, cervical stenosis, cervical/vaginal infection, uterine bleeding or menses, known pregnancy, invasive carcinoma of the cervix, recent uterine perforation, medical contraindication or intolerance to anesthesia.

Warnings: Failure to follow any instructions or to heed any Warnings or Precautions could result in serious patient injury. The HERizon Hysteroscopy System should only be used by physicians trained in hysteroscopy and hysteroscopic surgery using powered instruments. Healthy tissue can be injured, e.g., perforation. Use every available means to avoid such injury.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. The physician using the system must be trained in hysteroscopy. Prior to use, consult the Instructions for Use and/or Operators Manual. For Indications and a complete list of contraindications, potential adverse events, warnings and precautions.



MANUFACTURED BY
Jiangsu Jiyuan Medical Technology Co., Ltd.
NO.48 CECEP (Taizhou) Environmental Protection Technology Industry
Park, No.59 Meilan Road, Hailing District, Taizhou, Jiangsu, China.

www.minervasurgical.com
© 2024 Minerva Surgical and its affiliates.
All rights reserved

minerva®
The Uterine Health Company

K0309 Rev. B